


EZ Link User Guide

If you need any further assistance call the EZ Link Support Line 0800 439 5465
(0800 4 EZLINK)





Contents

Getting to know the Zebra Device

Power Button	3
Removing from Charger	3
Unlocking the device & Scanning Flasks	3
Home screen	3

EZ Link App

Logging in	4
Home page	4
Test setup	4

Herringbone

Record - Alternative Scan ON	5
Record - Alternative Scan OFF	5
Reverse record/scan	5
Add missed animal (Alternative Scan OFF)	6
Add missed flask (Alternate scan OFF)	6
Short row - skip all remaining empty bails	6
Clear flasks	6
Clear Row (Animal ID & Flask)	6

Rotary

Record	7
Reverse record/scan	7
Slow milking animal	7

Herringbone & Rotary

Enter animals	8
Enter animal ID via EID wand - supplied by LIC	8
Edit animal ID	8
Animals unknown to MINDA	8
Duplicate animals	9
OAD animals	9
Sick Animals	9
Adding/Editing an abnormal code	10
Re-scan flask	10
Duplicate flask	10
Overflow flask	11
Clear Bail (Delete animal and flask)	11
Query flask	11

Review

Review	12
--------	----

Getting to know your Zebra Device

Power button


Press and hold down once.

If the screen does not light up. Hold the Power Button down until the screen turns on.



Scanning

There are scanning buttons on either side of the device. There will also be a scan button on screen when required.



Scan →

Scanning lens

← Scan



Unlock Screen


Unlock device by swiping the screen upwards.

Remove from Charger

To unclip from the charger push one side of the charger clip outwards.

Home screen

If you find yourself at this screen, Tap on the blue EZ Link App symbol to take you back into the app.




EZ Link App

Logging in

1. Enter your participant code.

If you are unsure of your code refer to the sticker on device or EZ Link box.


2. Select 'Sign In'.



The login screen features the LIC logo at the top, followed by the text 'EZ Link'. Below this is a white input field labeled 'Participant code'. At the bottom is a dark blue button with the text 'Sign In'.

Home page


- Green tick indicates an action is complete
- If a field is greyed out (e.g. Record AM) it indicates the test setup is incomplete or the previous session does not have any animals recorded in it.
- Shows how many animals are recorded in each session.
- Review
 - o Look up and edit sampled animals.
- Query Flask
 - o Scan flask to find out which animal it is allocated to.
- Logout
 - o At end of session.



The home page shows the LIC logo and 'EZ Link' at the top. Below is a header with 'LIC Innovation Farm' and 'PDGY'. A list of actions follows: 'Test Setup' with a green checkmark, 'Record PM' with '7 animals' and a green checkmark, 'Record AM' with a blue arrow, and 'Review' with a blue arrow. At the bottom are 'Query Flask' and 'Logout' options.

Test setup





1. Fill out mandatory fields
 - Previous cups on date
 - Previous cups on time.
2. Optional settings
 - Enable EID wand - Yes/No (EID wands booked through LIC).
 - Alternate scan - On/Off (herringbone only)
 - o On - Enter animal ID, scan flask, enter animal ID, scan flask
 - o Off - Enter all animal IDs in the row, then scan all flasks in the row.
3. Select **NEXT** to save your changes.




The test setup screen has a back arrow and 'Test Setup' title. It shows 'LIC Innovation Farm' and 'PDGY'. The 'Test date' is set to 'Tue, 06 Jul 2021'. Below are two fields for 'Previous cups on date' and 'Previous cups on time', both with red borders and warning icons. At the bottom is an 'Enable EID Wand' dropdown set to 'No'.

Record - Alternative Scan ON






Enter animal ID, scan flask, enter animal ID, scan flask


1. Enter animal ID.
2. Scan flask.
3. Blue flask  indicates successful scan.
4. Continue entering animal IDs and scanning flasks until end of row.
5. Select  ADD EXTRA BAIL if there is an extra animal to record in the row.
6. To finish row scroll to bottom of record screen and select 
7. When finished session select  to go back to home screen.
8. Select 'Logout'.



Record - Alternative Scan OFF

Enter all animal IDs in the row, then scan all flasks in the row

1. Enter all animal IDs in the row.
2. Select  ADD EXTRA BAIL if there is an extra animal to record in the row.
3. Select 
4. Scan all flasks in row.
5. Blue flask  indicates successful scan.
6. To finish a row, scroll to bottom of record screen and select 
7. When finished session select  to go back to home screen.
8. Select 'Logout'



Reverse record/scan

Select  to change bail numbers from ascending to descending

If	Then
Alternative scan ON	Select before entering first animal ID in row
Alternative scan OFF	Select before entering any animal IDs or before scanning any flasks


Herringbone

Add missed animal (Alternative Scan OFF)

A missed animal can be added while recording animal IDs in a herringbone shed.



1. Tap on the bail where the missed animal should be.
2. Tap 'Add missed animal'.

This action moves all animals entered down one bail.



Add missed flask (Alternate scan OFF)

1. Tap on the 3 dots in top right hand corner of record screen.
2. Tap 'Add Missed Flask'.
3. Scan back through flasks until missing flask is located.
4. When missing flask is found users will be prompted to scan remaining flasks in the row.



Warning

Skip the remaining 14 bails and begin scanning?

CANCEL

BEGIN SCANNING

Short row - skip all remaining empty bails


1. Scroll to bottom of record screen.
2. Select **SCAN FLASKS**
3. Confirm by selecting 'Begin scanning' on warning.

Clear flasks

1. Tap on the 3 dots in top right hand corner of record screen.
2. Select 'Clear Flasks'.
3. All flasks scanned in row will be cleared.

Clear Row (Animal ID & Flask)

4. Tap on the 3 dots in the top right hand corner of record screen.
5. Select 'Clear Row'.
6. All animal IDs and flasks in row will be cleared.











Rotary

Record


Enter animal ID, scan flask, enter animal ID, scan flask

1. Enter animal ID.
2. Scan flask.
3. Blue flask  indicates successful scan.
4. Continue entering animal IDs and scanning flasks until the end of the session.
5. If there is an empty bail, select 
6. When finished session select  to go back to home screen.
7. Select 'Logout'.

 Record PM
 


1	10	
2	11	
3	<input type="text"/>	

 TAG #
 BIRTH ID


 OAD
  SICK

1	2	3
4	5	6
7	8	9
DELETE	0	ENTER

Reverse record/scan

Select  to change bail numbers from ascending to descending – Must be selected before any animals are recorded in the session.

Slow milking animal



If	Then
Animal goes around rotary with no change of flask	<ul style="list-style-type: none"> • Skip the bail or • Record the animal and flask again – EZ Link will recognise the animal as a slow milker
Animal goes around rotary and flask is changed	<ul style="list-style-type: none"> • Duplicate animal screen will appear with option to record second flask as overflow <div>  Mark as overflow </div>

Herringbone & Rotary

Enter animals

- Enter animal ID using keypad on device.
- Press **ENTER** to move to the next action.
- To record a Birth ID tap the button circled in red to the right.
- The format for recording Birth IDs is PTPT-YY-XXX.
 - o Example: PDGY-17-123
 - o The dash in between fields automatically generates - users do not need to add it.
- The participant field is pre-filled with the herd's participant code, to remove this use the **DELETE** button.

Enter animal ID via EID wand - supplied by LIC



1. Enable EID wand in Test Setup by changing from No to Yes.
2. Select **NEXT**
3. Go back into either 'Record AM' or 'Record PM'.
4. Tap  in top right hand corner of record screen to enable EID wand.
5. Icon will change when EID wand is connected 

Edit animal ID

1. Tap on the bail with the animal ID you wish to edit.
2. Select Edit animal ID.
3. Use the **DELETE** button on the keyboard to remove the incorrect number.
4. Re-type correct animal ID.
5. Select SAVE.



Animals unknown to MINDA®


- Herd details are pre-downloaded into EZ Link from MINDA.
- Animals not recorded in MINDA will show with  symbol.
- Duplicate animals designated with animal IDs of 10,001 or higher will also display with  symbol.
- Animals unknown to MINDA can be viewed and edited in the review screen.


Duplicate animals

If a duplicate animal is recorded the screen to the right will display.

- Select 'EDIT' to change an animal ID.
- If you mark current as correct the incorrect animal will display with 10,000 in front of the ID.

For example: If there is a duplicate of cow 50, the incorrect animal will display as 10050.

- Animals marked as duplicate will show as 'Unknown to MINDA' 
- If the duplicate is in the same bail users will be given the option to record an overflow flask.
- Duplicate animals can be viewed and edited in the review screen.



OAD animals

There are two ways to record OAD animals.

1. Tap on the OAD button once to mark one animal as OAD.
2. Tap on the OAD button twice to record multiple animals as OAD.

The icon will display with lock symbol  OAD


3. Confirm 'LOCK OAD'.
4. To turn off OAD, tap the button again.

Or

1. Tap on animal ID to bring up edit options.
2. Select 'Mark as OAD'.

To **unmark** as OAD

1. Tap on animal ID to bring up edit options.
2. Select 'Unmark as OAD'.




Confirm Lock OAD

Animals from this point will be marked as OAD until unlocked

CANCEL

LOCK OAD



Sick Animals

There are two ways to record an animal as sick.

1. Tap on the SICK button once to mark one animal as sick.
2. Tap on the SICK button twice to record multiple animals as sick.

The icon will display with lock symbol  SICK

3. Confirm 'LOCK SICK' on.
4. To turn off SICK, tap again.

Or follow the steps for Adding/Editing an abnormal code on page 10.



Confirm Lock Sick

Animals from this point will be marked as SICK until unlocked

CANCEL


LOCK SICK



Herringbone & Rotary


Adding/Editing an abnormal code

1. Tap on animal ID to bring up edit options
2. Select 'Edit abnormal code'
3. Select appropriate abnormal code
4. Select 'DONE'










Re-scan flask

1. Tap on animal ID to bring up edit options.
2. Select 'Re-scan flask'.




Duplicate flask


If	Then
The duplicate flask is in the same row (herringbone) or rotation (rotary).	Users are given the below options.  Rescan previous flask  Rescan current flask  Mark samples as mixed  Undo last action
The duplicate flask is NOT in the same row (herringbone) or rotation (rotary).	Users are given the below options.  Rescan current flask  Mark samples as mixed  Undo last action

Herringbone & Rotary

Overflow flask


There are 3 ways users can record an overflow flask.

1. Tap on animal ID to bring up edit options - select Add overflow flask.
2. In a rotary if an animal goes around twice a pop-up will show she is a duplicate and users will be given the option to mark second sample as overflow 
3. Review screen - tap on animal ID - users will be given the option to add overflow flask for either the PM or AM session.



Clear Bail (Delete animal and flask)


1. Tap on animal ID to bring up edit options - select clear bail.
2. User will be asked to confirm.
3. Animal ID and flask will be deleted.



Query flask

Use when uncertain whether a flask has been scanned.


1. Find query flask on home screen or tap on 3 dots in top right hand corner of record screen.
2. Select 'Query Flask'.
3. Scan the flask in question.
4. Animal ID of linked flask will display - or will display 'no record found'.
5. User is then able to scan another flask (if required).





Review

Users are now able to review throughout a milking session

- Search by entering animal ID next to 
- Edit by tapping on the animal ID.
- Users can filter by the following:
 - o Duplicate animals
 - o Unknown to MINDA
 - o No calving date
 - o Single-sample animals
 - o Lactating animals not sampled
 - o OAD
 - o Abnormal code animals.

← Review		
<input type="text" value="Q"/>		
Filter by ▾		
	PM 3	AM 3
1		
3		
4		
10		

If you need any further assistance call the EZ Link Support Line 0800 439 5465 (0800 4 EZLINK)